

Using the New ISBT 128 Product Lookup Program Website Application: Exporting Search Results

The search results can be exported from the ISBT 128 Product Lookup Program Website Application which can be imported into another program.

Exporting All Inclusive Results:

1. Click the “Export All Inclusive Results” button. This downloads all search results into a Microsoft Excel Comma Separated Values (CSV) File.
2. Locate the downloaded file. Right-click the file and select open with “Notepad”.
3. Click “File” and “Save As”. Save the results as a Text Document (.txt).
 - This semicolon-delimited text file can be used to import into another program.

See the following images.

Search by Product Description

Category

Cellular Therapy

Subcategory

Cellular Therapy

Class

HPC, APHERESIS

Attribute Group

Mobilization

Attribute Value

Mobilized

Selected Attributes

Citrate/XX/<=-150C
V0002035
5% DMSO
V0045004
Cryopreserved
V0050005
Mobilized
V0051002



Clear Attributes

1

Export All Inclusive Results

Clear ALL values

Using the New ISBT 128 Product Lookup Program Website Application: Exporting Search Results




Exporting Selected Results:

1. Scroll down in the program to find the inclusive results listed.
2. Select the preferred Product Description Codes to export.
3. Click the "Export" button after making selections. This downloads all selected search results into a Microsoft Excel Comma Separated Values (CSV) File.
4. Locate the downloaded file. Right-click the file and select open with "Notepad".
5. Click "File" and "Save As". Save the results as a Text Document (.txt).
 - This semicolon-delimited text file can be used to import into another program.

See the following images.

Using the New ISBT 128 Product Lookup Program Website Application: Exporting Search Results

All Inclusive Results  1

Selected: 3

☒ Select all on this page











☐ Un-select this page

 Export

 Copy

Filter



3

Export	Code	Description	Actions
<input checked="" type="checkbox"/>	S1556	HPC, APHERESIS Citrate/XX/<=-150C 5% DMSO 3rd Party Comp:Yes Other Additives:Yes Cryopreserved Mobilized	
<input type="checkbox"/>	S3114	HPC, APHERESIS Citrate/XX/<=-150C 5% DMSO Concurrent plasma + other Cryopreserved Mobilized	
<input checked="" type="checkbox"/>	S3299	HPC, APHERESIS Citrate/XX/<=-150C 5% DMSO 3rd Party Comp:Yes Other Additives:Yes Cryopreserved Mobilized Plasma reduced	
<input checked="" type="checkbox"/>	S3401	HPC, APHERESIS Citrate/XX/<=-150C 5% DMSO Cryopreserved Mobilized	
<input type="checkbox"/>	S3404	HPC, APHERESIS Citrate/XX/<=-150C 5% DMSO Other Additives:Yes Cryopreserved Mobilized Plasma reduced	
<input type="checkbox"/>	S3477	HPC, APHERESIS Citrate/XX/<=-150C 5% DMSO Concurrent plasma Cryopreserved Mobilized	
<input type="checkbox"/>	S3478	HPC, APHERESIS Citrate/XX/<=-150C 5% DMSO Other Additives:Yes Cryopreserved Mobilized	
<input type="checkbox"/>	S3483	HPC, APHERESIS Citrate/XX/<=-150C 5% DMSO 3rd Party Comp:Yes Cryopreserved Mobilized	
<input type="checkbox"/>	S3484	HPC, APHERESIS Citrate/XX/<=-150C 5% DMSO 3rd Party Comp:Yes Concurrent plasma Cryopreserved Mobilized	
<input type="checkbox"/>	S3494	HPC, APHERESIS Citrate/XX/<=-150C 5% DMSO Cryopreserved Mobilized Plasma reduced	

Show 50 of 10 rows per page

Page 1 of 1

Using the New ISBT 128 Product Lookup Program Website Application: Exporting Search Results



For any questions about using the new ISBT 128 Product Lookup Program Website Application, please contact the ICCBBA office:

- iccbba@iccbba.org
- +1 909 793 6516